

ANANTH


Anjiyografik Diagnostik Kateteri

Sipariş Bilgisi

S.No.	Uzunluk	French Boyu	Kıvrım Şekli	Giriş Yolu	Kıvrım Ayrıntısı	Ürün Kodu
1	100 cm / 110 cm	5 Fr	Tiger Straight	Radial	3.50 cm	ATS- 3500 – 5Fr
2					4.00 cm	ATS- 4000 – 5Fr
3					4.50 cm	ATS- 4500 – 5Fr
4			Tiger Angled	Radial	3.50 cm	ATA- 3500 – 5Fr
5					4.00 cm	ATA- 4000 – 5Fr
6					4.50 cm	ATA- 4500 – 5Fr
7			Judkins Left	Femoral	3.50 cm	AJL- 3500 – 5Fr
8					4.00 cm	AJL- 4000 – 5Fr
9					4.50 cm	AJL- 4500 – 5Fr
10			Judkins Right	Femoral	3.50 cm	AJR- 3500 – 5Fr
11					4.00 cm	AJR- 4000 – 5Fr
12					4.50 cm	AJR- 4500 – 5Fr
13			Pig Tail	Femoral	90°	APT- 0900 – 5Fr
14					145°	APT- 1450 – 5Fr
15					155°	APT- 1550 – 5Fr
16			Amplatz Left	Femoral	AL1	AAL- AL01 – 5Fr
17					AL2	AAL- AL02 – 5Fr
18					AL3	AAL- AL03 – 5Fr
19			Amplatz Right	Femoral	AR1	AAR- AR01 – 5Fr
20					AR2	AAR- AR02 – 5Fr
21					AR3	AAR- AR03 – 5Fr
22	100 cm / 110 cm	6 Fr	Tiger Straight	Radial	3.50 cm	ATS- 3500 – 6Fr
23					4.00 cm	ATS- 4000 – 6Fr
24					4.50 cm	ATS- 4500 – 6Fr
25			Tiger Angled	Radial	3.50 cm	ATA- 3500 – 6Fr
26					4.00 cm	ATA- 4000 – 6Fr
27					4.50 cm	ATA- 4500 – 6Fr
28			Judkins Left	Femoral	3.50 cm	AJL- 3500 – 6Fr
29					4.00 cm	AJL- 4000 – 6Fr
30					4.50 cm	AJL- 4500 – 6Fr
31			Judkins Right	Femoral	3.50 cm	AJR- 3500 – 6Fr
32					4.00 cm	AJR- 4000 – 6Fr
33					4.50 cm	AJR- 4500 – 6Fr
34			Pig Tail	Femoral	90°	APT- 0900 – 6Fr
35					145°	APT- 1450 – 6Fr
36					155°	APT- 1550 – 6Fr
37			Amplatz Left	Femoral	AL1	AAL- AL01 – 6Fr
38					AL2	AAL- AL02 – 6Fr
39					AL3	AAL- AL03 – 6Fr
40			Amplatz Right	Femoral	AR1	AAR- AR01 – 6Fr
41					AR2	AAR- AR02 – 6Fr
42					AR3	AAR- AR03 – 6Fr

ANANTH

Anjiyografik Diagnostik Kateteri



Relisys Medical Devices Limited

Survey.No.312, Pocharam Road, Mangalpally Village,
Ibrahimpatnam Mandal, Ranga Reddy District,
Telangana State - 501 510, INDIA
Email: info@relisysmedicaldevices.org
website: www.relisysmedicaldevices.org
Tel: +9140 23493349, Fax: +9140 23493301



Relisys Medical Devices Ltd.
An ISO 13485 Certified Company

ANANTH

Anjiografik Diagnostik Kateteri

- Sol ve sađ koroner arterin radyal ve femoral ulaşımı için **Hidrofil kaplı** 5 Fr ve 6 Fr kateter tasarımı
- Seçkin kıvrım şekli hafızası
- Üstün torklanabilirlik
- Çok iyi kıvrım direnci
- Mükemmel radyoopak izlenebilirliği
- Damar hasarını azaltmak için atravmatik yumuşak uç
- Yüzey bitişinde arteriyal spazmları azaltmak için son teknoloji
- Kateter 950 psi'a kadar yüksek basınca direnebilir
- 0.038" kılavuz tel uyumluluđu

